

03

# FUNGI KIT

*Plotting Guide*

## Step 1: Saving the Data from the Raspberry Pi Pico

! Disconnect the **analogue-to-digital converter (ADC)** from the Raspberry Pi Pico by unplugging the wires before continuing.

⇒ Reopen Thonny, then plug the Raspberry Pi Pico back into the computer and press the red **Stop** button to tell **Thonny** to find it.





fig. 1

⇒ We want to save the data file from the Raspberry Pi Pico to a simple location. In the **Files** tab on the left of **Thonny**, click on the **C:\** (fig. 2).




fig. 2

⇒ In the Files window, navigate to **Users** → **[Your Computer Name]** → **Documents** → **FungiKit**. You should now see a file named **fungi\_kit\_plotter.py**.

⇒ **Note:** You won't see the text "**Your Computer Username**" — it will show your actual computer name instead.

- Right-click on the **data.csv** file and click “**Download to C:\Users\[Your Computer Name]\Documents\FungiKit**”.
- Also open the file **fungi\_kit\_plotter.py** file by double clicking it.




fig. 3

## Step 2: Configure the Interpreter

- Go to the top toolbar and select **Run > Configure Interpreter...**




fig. 4

⇒ In the first drop-down, select **Local Python 3** and then click **OK**.




fig. 5

### Step 3: Plotting the Data

⇒ In the **fungi\_kit\_plotter.py** file you opened earlier, scroll down and rename the **"TITLE"**. Make sure to include the date of the recording taken in your title.

```
28
29 # Plot each of the sensor readings (A0-A3) over time
30 plt.plot(data["timestamp"], data["A0"], label="A0")
31 plt.plot(data["timestamp"], data["A1"], label="A1")
32 plt.plot(data["timestamp"], data["A2"], label="A2")
33 plt.plot(data["timestamp"], data["A3"], label="A3")
34
35 # Label the x-axis (horizontal) and y-axis (vertical)
36 plt.xlabel("Time")
37 plt.ylabel("Value (V)")
38
39 # Give the graph a title
40 plt.title("TITLE")
41
42 # Show a small box that explains which line is which
43 plt.legend()
44
45 # Display the graph on the screen
46 plt.show()
```

fig. 6

⇒ Click the **Green Play Button** in the top left corner of **Thonny** to run the program.




fig. 7

⇒ After a few seconds, a new window will open with your data plot



fig. 8

⇒ Congratulations! You have now finished recording and plotting your data. Your teacher and the tutors will now explain how to analyse the recording.

⇒ If you have finished this task quickly, flip over to the next page for an extension task.

## Extension Task: Smoothing the Data with Filtering

- ☞ You may notice that the data you plotted is quite **noisy**, which means it jumps up and down a lot instead of following a smooth trend. Noisy data can happen for many reasons, like small measurement errors or natural variations in the mushrooms. To better see the overall pattern, we can use a method called **filtering** to smooth the data. Filtering reduces the random ups and downs, making it easier to spot trends and understand what the data is really showing.
- ☞ To filter the data, we need to add a few more lines of code **before plotting the data**. Add the following highlighted lines of code (fig. 9).

```

26 # Drop (remove) any rows that have missing values in these columns
27 data = data.dropna(subset=["timestamp", "A0", "A1", "A2", "A3"])
28
29 for col in ["A0", "A1", "A2", "A3"]:
30     data[col] = data[col].rolling(window=30, center=True).mean()
31
32 # Plot each of the sensor readings (A0-A3) over time
33 plt.plot(data["timestamp"], data["A0"], label="A0")
34 plt.plot(data["timestamp"], data["A1"], label="A1")
35 plt.plot(data["timestamp"], data["A2"], label="A2")
36 plt.plot(data["timestamp"], data["A3"], label="A3")
37
38 # Label the x-axis (horizontal) and y-axis (vertical)
39 plt.xlabel("Time")
40 plt.ylabel("Value (V)")

```

fig. 9

- ☞ Run the script again using the **Green Play Button**.
- ☞ You should see that the data on the new plot is smoothed, and the spiking activity is more clearly visible.




fig. 10

- ☞ Experiment with different window sizes by changing the number in **[window=VALUE]** and observe how it changes the graph.